



INTRODUCTION

Next Generation Sequencing (NGS) library construction is a workflow bottleneck. The processes of library construction for DNA and RNA sequencing are complex, error prone, costly, and in general are time consuming. Current commercial products offer multistep, multi-hour sample processing that can be cost and labor prohibitive. These issues have created a demand for a simple, rapid, and cost-effective library construction products that offer users application flexibility while minimizing construction complexity for low to ultra-high throughput sample processing. SeqOnce has developed a novel library technology that is rapid and minimizes sample processing complexity. The five tube kit is comprised of formatted master mixes for a simple and stable user workflow. The 12 minute library construction uses a single master-mix that when combined with fragmentation and PCR steps produces libraries in less than 50 minutes. A single size selection step occurs after PCR. The PCR free and pre-fragmented DNA protocols are less than 30 minutes.

The kit contains three master mixes and includes adapters, which maximally promotes product flexibility across multiple applications. With master-mixes that are stable for multiple days at ambient temperatures, automation on liquid handling platforms is effortless due to the simple workflow. The technology has been validated on a variety of human sample types and applications with input ranges of 10ng- 100ng. The sequencing data is either equivalent or superior to other competing library products, with high mapping rates and excellent performance across variable GC content. The technology is currently optimized for the Illumina platform and alternative NGS platforms are under evaluation for future development.

WORKFLOW & COMPONENTS



Single Tube Additive Reagent Library Prep

- Libraries in under 50 minutes
- Two step protocol, three steps with PCR
- Reagents stable for >24 hours at ambient temperature
- Minimal setup due to master mix format
 - Five reagent tube kit:
 - Three master mix tubes for:
 - Fragmentation
 - Library Construction
 - PCR
 - Two remaining tubes:
 - Stock Adapter
 - Adapter Dilution Buffer



RHINOSEQ REAGENTS

Fragmentation Master Mix Characteristics:

- Single volume for 5 – 100ng inputs
- Adjustable for controlled fragment distribution
- No buffer/enzyme pipetting
- Single tube buffer/enzyme mix
- Identical conditions for microbial human gDNA fragmentation

FRAGMENTATION MASTER MIX REAGENT SINGLE TUBE – SINGLE VOLUME

Fragmentation Distribution




Figure 5: Histogram profile of Rhinoseq short or long fragmentation of Human DNA. Various input amounts of Human gDNA were fragmented using the Rhinoseq kit to generate short (left) or long (right) fragment distributions. After fragmentation, samples were 2X SPRI[®] bead cleanup and analyzed using an Agilent High Sensitivity DNA Assay[™].

Library Fragment Length Distribution




Figure 6: Histogram profile of Rhinoseq short or long insert Human gDNA Libraries. Nucleotide content over a 150bp window of sequencing read illustrates similar start-site complexity between Covaris (left) and random enzymatic fragmentation (right) using the Rhinoseq kit.

LIBRARY PREP MASTER MIX REAGENT SINGLE TUBE – SINGLE VOLUME – 12 MIN PREP

Minimal Start Site/Coverage Bias




Figure 7: Read start-site bias for Rhinoseq E. coli Libraries. Nucleotide content over a 150bp window of sequencing read illustrates similar start-site complexity between Covaris (left) and random enzymatic fragmentation (right) using the Rhinoseq kit.

MICROBIAL SEQUENCING

WHOLE GENOME Third Party- Independent Tester

Genome Coverage & GC Bias

Species	GC %	Library Input	Coverage Depth	Genome Coverage (unique reads)
<i>Staphylococcus</i>	32.7	30	10 - 20	+93%
<i>Escherichia coli</i>	51	10, 20, 30, 40	19 - 24	+96%
<i>Rhodobacter</i>	70	30	9.5 - 32	+96%

Table 1: Bacterial sequencing and coverage statistics. Genome structure and sequencing statistics of *Staphylococcus*, *E. coli*, and *Rhodobacter*.




Figure 8: GC Bias Comparison for Microbial Samples. Rhinoseq (grey), *E. coli* (blue), and *Staphylococcus* (yellow) GC bias was assessed by using Picard CollectGCBiasMetrics. A normalized coverage value of 1 is representative of an absence of sequencing bias. The Rhinoseq kit results in minimal GC bias in AT or GC rich genomes.

Community Characteristics

Species	GC %	Gram Strain	gDNA Abun. (B)
<i>Pseudomonas aeruginosa</i>	66.2	-	12
<i>Escherichia coli</i>	56.8	-	12
<i>Salmonella enterica</i>	52.2	-	12
<i>Lactobacillus fermentum</i>	52.8	+	12
<i>Enterococcus faecalis</i>	37.5	+	12
<i>Staphylococcus aureus</i>	32.7	+	12
<i>Listeria monocytogenes</i>	38.0	+	12
<i>Bacillus subtilis</i>	43.8	+	12
<i>Saccharomyces cerevisiae</i>	38.4	Yeast	2
<i>Cryptococcus neoformans</i>	48.2	Yeast	2

Table 2: Community Standard Characteristics The Zymo Research Microbial Community DNA Standard was used to assay library construction efficiency using the SeqOnce Rhinoseq kit. The standard contains a defined and characterized composition of 5 Gram-Positive and 3 Gram-Negative bacteria plus 2 yeast species. The GC content of the standard range is 15%-85%.

Community Sequencing




Figure 9: Taxonomic Identification The Zymo Research Microbial Community DNA Standard was used to assay library construction efficiency using the Rhinoseq kit. Profiling was with 2x75bp shotgun sequencing (1.5-2 million reads, Zymo: 176 million reads). Left panel, data analysis was conducted by reference read mapping and counting. Right panel, microbial composition was determined by using mOTU (<http://www.bork.embl.de/software/motu/>). Relative abundance deviation in average: <15% (indicated by white not included in the analysis).

HUMAN WGS

DEEP COVERAGE Joseph Boland – Independent Tester NCI/NIH, Leidos Biomedical Research Inc.




Figure 10: Sequence Quality Base phred quality (top) and mean read quality scores (bottom) show high quality data from 50ng and 100ng NA12878 and NA12892 DNA libraries.

Figure 11: Human Chr. 1 Coverage. Coverage across Chr. 1 for NA12878 and NA12892 50ng (blue) and 100ng (pink) input libraries.

Figure 12: ROC Plots. True positive calls (y-axis) versus false positive calls (x-axis) for 50ng and 100ng NA12878 and NA12892 DNA libraries are plotted. Fifty nanogram libraries are in blue and 100ng libraries are in orange.

SHALLOW COVERAGE




Figure 13: Mean Coverage Per Chromosome. The mean coverage per chromosome was consistent, ranging between 3-4x coverage.

GC Bias




Figure 14: Rhinoseq Minimal GC Bias The Rhinoseq kit resulted in minimal GC bias despite sequencing to shallow depth. GC bias was assessed by using Picard CollectGCBiasMetrics.

HUMAN FFPE




Figure 15: Mild FFPE Library Pre- and post- library prep was assayed. The pre-fragmented DNA was constructed with the standard Rhinoseq kit.

Figure 16: Severe FFPE Library Pre- and post- library construction was assayed. The pre-fragmented DNA was constructed with the standard Rhinoseq kit and produced libraries in 30m.

Mapping

	GENOME TERRITORY	HUMAN COV	GENOME COV	GENOME COV
Severe-1	3066216270	0.703822	1.508231	90.40%
Severe-2	1176392479	1.834482	1.963387	
Severe-3	3066216270	0.55336	1.315136	90.40%
NON-ZERO	2754760736	18.403834	8.187294	

HLA TYPING

Independent Tester – HLA Typing Company



Figure 17: HLA Amplicon Libraries Plot represents a subset of Bioanalyzer traces of 48 libraries generated with the Rhinoseq kit using 22-53ng HLA multiplexed amplicons (4-7kb).

Figure 18: HLA Coverage Per Amplicon Each amplicon had 100% coverage and 99.7% concordance.

FUTURE PRODUCTS



CONCLUSION

The RhinoSeq kit allows for DNA library prep that is both straightforward and rapid. Total library processing time (less than 50 minutes) is faster than other commercial kits (three to four times faster). The simplicity of the five tube two-step protocol, three-step with PCR, enables experienced and inexperienced users to easily prepare libraries for sequencing. The flexibility of this system also allows users to readily scale the number of samples processed, easy to automate, and can be applied to a variety of sequencing applications.

1. Bioinformatics was performed by independent testers whose identities are indicated on the Bioinformatics page.
2. Rhinoseq kit was used for sequencing of human DNA. The sequencing data was analyzed using Picard CollectGCBiasMetrics. The GC content of the standard range is 15%-85%.
3. The Genome Analysis Toolkit (GATK) is a software package for analyzing next-generation DNA sequencing data. Version 4.0.6. © 2012 GENOME RESEARCH INTERNATIONAL. ALL RIGHTS RESERVED.
4. Picard is a software package for analyzing next-generation DNA sequencing data. Version 2.19.0. © 2012 GENOME RESEARCH INTERNATIONAL. ALL RIGHTS RESERVED.
5. The Zymo Research Microbial Community DNA Standard was used to assay library construction efficiency using the SeqOnce Rhinoseq kit. The standard contains a defined and characterized composition of 5 Gram-Positive and 3 Gram-Negative bacteria plus 2 yeast species. The GC content of the standard range is 15%-85%.
6. mOTU is a software package for analyzing next-generation DNA sequencing data. Version 1.0. © 2012 GENOME RESEARCH INTERNATIONAL. ALL RIGHTS RESERVED.
7. Picard is a software package for analyzing next-generation DNA sequencing data. Version 2.19.0. © 2012 GENOME RESEARCH INTERNATIONAL. ALL RIGHTS RESERVED.