

AzureSeq One-Step Universal RT-qPCR Kit SARS-CoV-2

Instructions for Use

Intended Use:

The AzureSeq One-Step Universal RT-qPCR Kit SARS-CoV-2 is an RT-qPCR test intended for the qualitative detection of nucleic acid from the 2019-nCoV in nasopharyngeal (NP) and oropharyngeal (OP) swabs from individuals with signs and symptoms of infection who are suspected of COVID-19. Testing is limited to laboratories - certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, to perform high complexity tests, or by similarly qualified non-U.S. laboratories.

Results are for the identification of 2019-nCoV RNA. The 2019-nCoV RNA is generally detectable in nasopharyngeal and oropharyngeal swabs during the acute phase of infection. Positive results are indicative of active infection. Laboratories within the United States and its territories are required to report all positive results to the appropriate public health authorities.

Negative results do not preclude 2019-nCoV infection and should not be used as the sole basis for patient management decisions. Negative results must be combined with clinical observations, patient history, and epidemiological information.

The AzureSeq One-Step Universal RT-qPCR Kit SARS-CoV-2 is intended for use by qualified, trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and in vitro diagnostic procedures. The AzureSeq One-Step Universal RT-qPCR Kit SARS-CoV-2 is only for use under the Food and Drug Administration's Emergency Use Authorization.

Test Principle:

Nucleic acids are isolated and purified from nasopharyngeal and oropharyngeal swabs using a previously FDA cleared nucleic acid extraction system. Sample input and elution volumes are system dependent. The purified nucleic acid is reverse transcribed into cDNA by combining the nucleic acid with the AzureSeq One-Step Universal RT-qPCR Kit SARS-CoV-2 master mix. In the process, the probe anneals to a specific target sequence located between the forward and reverse primers. During the extension phase of the PCR cycle, the 5' nuclease activity of Taq polymerase degrades the probe, causing the reporter dye to separate from the quencher dye, generating a fluorescent signal. With each cycle, additional reporter dye molecules are cleaved from their respective probes, increasing the fluorescence intensity. Fluorescence intensity is monitored at each PCR cycle.



Materials Provide and Storage:

Concentration: 4X

Storage and Handling: Store at -20°C upon arrival

Ordering Information:

| Item Number | Total Volume Received | Quantity Received | Total number of reactions which can be obtained when using the following reaction sizes | |
|---------------|---|--|---|---|
| | | | 20µL Reactions | 10µL Reactions |
| QSP-KIT-200 | ITP-MM-200-PC: 1mL, RT-200-PC: 100uL, DTT-100-200uL-PC: 200uL | ITP-MM-200-PC: 1x1mL, RT-200-PC: 1x100uL, DTT-100-200uL-PC: 1x200uL | 200 | 400 |
| QSP-KIT-1000 | ITP-MM-200-PC: 5mL, RT-200-PC: 500uL, DTT-100-200uL-PC: 1000uL | ITP-MM-200-PC: 5x1mL, RT-200-PC: 5x100uL, DTT-100-200uL-PC: 5x200uL | 1000 | 2000 |
| QSP-KIT-2000 | ITP-MM-200-PC: 10mL, RT-200-PC: 1000uL, DTT-100-200uL-PC: 2000uL | ITP-MM-200-PC: 10x1mL, RT-200-PC: 10x100uL, DTT-100-200uL-PC: 10x200uL | 2000 | 4000 |
| QS-SR-100*** | 5X RTScript Buffer-1x0.4mL, 10mM dNTP Mix-1x0.1mL, 100mM DTT Solution-1x0.1mL, 100µM Oligo-(dT)20 primer- 1x0.05mL, 100µM Random Hexamers-1x0.05mL, 40 Units/uL Rnase Inhibitor-1x0.01mL, 10ng/uL Positive Control RNA-1x0.01mL, Rnase-free Water-1x1.2mL | | 100 | 200 |
| EVA-300uL**** | EvaGreen Dye, 20X in Water; EVA: 2x150uL | | 200 | 400 |
| ROX-300uL**** | ROX Reference Dye, 20X in Water; ROX: 2x150uL | | Determined my volume required in assay by machine | Determined my volume required in assay by machine |

***This item can be purchased separately to perform a Two-Step Protocol: First Strand cDNA Synthesis followed by Amplification.

****This item can be purchased separately if dyes are required for qPCR assays.

Workflow:

Protocol for One-Step: The following reaction set up and general cycling conditions are recommended but can vary depending on the template and primers being used.

Reaction set-up for 20µL reaction volume Table 1:

| Component | Volume | Final Concentration |
|---|-------------|---------------------|
| 4X InhibiTaq Plus HotStart qPCR MasterMix | 5 µl | 1X |
| RTScript™, 200U/µL [#] | 0.5 µl | 5 units/µl |
| Target Specific Primers and Probe | 1µl | 1X |
| 100mM DTT [#] | 1µl | 5mM |
| Sample | X µl | 100 ng to 1 pg |
| EvaGreen Dye, 20X in Water (if required) | 1 µl | 1X |
| ROX reference dye, 25µM (See Table 2) | See Table 2 | See Table 2 |

| | | |
|-------------------------------|------------|------|
| Nuclease Free Water to volume | X μ l | N.A. |
| Total | 20 μ l | |

#Reactions need to be placed in cycler immediately following completion to maintain the integrity of the DTT and the RTScript™.

Table 2: ROX concentration recommendation for different instrument

| Type | Instrument | | ROX final concentration |
|----------|------------|--|-------------------------|
| | Company | Instrument Name | |
| No ROX | Roche | LightCycler 480, LightCycler 2.0 | None |
| | BioRad | iCycler, MyiQ, MiQ 2, iQ 5, CFX-96, CFX-384, Chromo4, MJ Opticon, Option2, MiniOpticon | |
| | Qiagen | Roto-Gene Q, Roto-Gene3000, Roto-Gene 6000 | |
| | Illumina | Eco RealTime PCR System | |
| | Eppendorf | Mastercycler realplex | |
| | Cepheid | SmartCyler | |
| Low ROX | ABI | 7500, 7500 Fast | 30nM |
| | Stratagene | MX4000P, MX3000P, MX3005P | |
| High ROX | ABI | 5700, 7000, 7300, 7700, 7900, 7900HT, 7900HT Fast, StepOne, StepOne plus, Viiia7 | 300nM |

Thermal cycling conditions:

Table 3: Recommended Cycling Conditions

| Cycling Step | Stage | No. of Cycles | Temperature | Holding Time |
|-------------------|-------|---------------|-------------|--------------|
| UNG Incubation | 1 | 1 | 25°C | 2 minutes |
| RT Incubation | 2 | 1 | 55°C | 15 minutes |
| Enzyme Activation | 3 | 1 | 95°C | 2 minutes |
| Amplification** | 4 | 40 | 95°C | 10 seconds |
| | | | 60°C | 60 seconds |

**Temperature and Holding Time will be based off primer set used.

Protocol for Two-Step: The following reaction set up and general cycling conditions are recommended but can vary depending on the template and primers being used.

Step 1A: cDNA synthesis without denaturation. General guidelines per 20 μ l reaction for setup without sample denaturation include, See Table 4

Table 4: cDNA synthesis without denaturation.

| Component | Stock Conc. | Final Conc. | 20uL Assay |
|------------------|-------------|---|--------------------------------|
| RNase-free water | - | - | Up to 20 μ L |
| RNA template | - | Total RNA: 10 pg-5 μ g or mRNA: 10pg-500ng | X μ L |
| Primer | 100 μ M | Gene specific primer: 10-20pg (50-100ng) Oligo-dT ₂₀ primer or random: 50pmol | 0.1-0.2 μ L 0.5 μ L |

| | | | |
|---------------------------------|----------------|------------|--------|
| 5X RTScript Buffer | 5X | 1X | 4µL |
| dNTP Mix | 10mM | 500µM each | 1 µL |
| DTT stock solution | 100mM | 5mM | 1 µL |
| RNase Inhibitor | 40 units / µL | 20 units | 0.5 µL |
| RTScript™ Reverse Transcriptase | 200 units / µL | 100 units | 0.5 µL |

Use reaction mix immediately to preserve the integrity of the enzyme.

Step 1B (optional): cDNA synthesis with denaturation.

- First, prepare the Template/Primer mix using the Table 5. Incubate the mix for 5 min at 65-70°C.

Table 5: Template/Primer Mix

| Component | Stock Conc. | Final Conc. | 20uL Assay |
|------------------|-------------|--|---------------------|
| RNase-free water | - | - | Up to 10µL |
| RNA template | - | Total RNA: 10 pg-5 µg or mRNA: 10pg-500ng | X µL |
| Primer | 100µM | Gene specific primer: 10-20pg (50-100ng) Oligo-dt ₁₅₋₂₅ primer or random: 50pmol | 0.1-0.2 µL 0.5µL |

- Second, prepare a Reaction mix using Table 6.

Table 6:

| Component | Stock Conc. | Final Conc. | 20uL Assay |
|---------------------------------|----------------|-------------|------------|
| RNase-free water | - | - | Up to 10µL |
| 5X RTScript Buffer | 5X | 1X | 4µL |
| dNTP Mix | 10mM | 500µM each | 1 µL |
| @DTT stock solution | 100mM | 5mM | 1 µL |
| %RNase Inhibitor | 40 units / µL | 20 units | 0.5 µL |
| RTScript™ Reverse Transcriptase | 200 units / µL | 100 units | 0.5 µL |

@Adding up to 5mM DTT may increase the yield and is recommended for individual optimization. %Addition of 20-40 Units of RNase inhibitor per assay is recommended when using low amounts of starting RNA. #100 Units of enzyme is recommended for standard assays, but increased transcription levels may be achieved with increasing the amount of enzyme up to 200 units. **Use reaction mix immediately to preserve the integrity of the enzyme.**

- Third, Add 10µl of Reaction Mix to 10uL of Template/primer mix and pipette gently up and down on ice.

Step 2: Incubation

Gene Specific primers: Incubate Reaction Mix for 30-60min at 50°C.

Oligo-dT or Random primers: Incubate Reaction Mix for 10min at 42°C followed by 30-60min at 50°C.

Step 3 (optional): Heat inactivation

Heat the mixture for 70°C for 10 min to inactivate Reverse Transcriptase.

Step 4(optional): RNA removal

Add 2 Units DNase-free RNase and incubate at 37°C for 20min.

The cDNA can now be used as a template in qPCR.

Reaction set-up for 20uL reaction volume

Table 7: Reaction Set Up

| Component | Volume | Final Concentration |
|---|-------------|---------------------|
| 4X InhibiTaq Plus HotStart qPCR | 5 μ l | 1X |
| Target Specific Primers and Probe | 1 μ l | 1X |
| 100mM DTT [#] | 1 μ l | 5mM |
| cDNA Template from Two-Step Protocol | X μ l | 100 ng to 1 pg |
| EvaGreen Dye, 20X in Water (if required) | 1 μ l | 1X |
| ROX reference dye, 25 μ M (See Table 2) | See Table 2 | See Table 2 |
| Nuclease Free Water to volume | X μ l | N.A. |
| Total | 20 μ l | |

[#]Reactions need to be placed in cyclers immediately following completion to maintain the integrity of the DTT.

Thermal cycling conditions:

Table 8: Recommended Cycling Conditions

| Cycling Step | Stage | No. of Cycles | Temperature | Holding Time |
|-------------------|-------|---------------|-------------|--------------|
| UNG Incubation | 1 | 1 | 25°C | 2 minutes |
| RT Incubation | 2 | 1 | 55°C | 15 minutes |
| Enzyme Activation | 3 | 1 | 95°C | 2 minutes |
| Amplification** | 4 | 40 | 95°C | 10 seconds |
| | | | 60°C | 60 seconds |

**Temperature and Holding Time will be based off primer set used.

Risk and Safety Information:

Safety Data Sheets are available online at www.seqonce.com

The user should carefully read all warnings, instructions or Safety Data Sheets provided by the supplier for any additional materials or chemicals required for the use of the AzureSeq One-Step Universal RT-qPCR Kit SARS-CoV-2. The user should also follow general safety precautions when handling biohazards, chemicals and other materials.

General Precautions:

1. Treat all samples, materials and instrumentation as potentially infectious.
2. Avoid potential contamination by employing good laboratory practices, wearing proper personal protective equipment, and decontaminating workspaces before and after use.
3. All instruments must be maintained and operated according to manufacturer's instructions.
4. Dispose of waste according to state and local regulations.
5. Only use DNase and RNase free consumables.
6. Do not use reagents beyond their expiration.

Manufacturing Information:










Empirical Bioscience
 2007 East Castle Drive, SE
 Grand Rapids, MI 49508

Support:

Phone: +1-626-460-3050
 Email: info@seqonce.com

Explanation of Symbols:

| Symbol | Definition |
|---|---|
|  | <i>In vitro</i> diagnostic medical device |
|  | Lot Number |
|  | Manufacturer |
|  | Date of Manufacture |
|  | Storage & Transport Conditions |

| | |
|---|-------------------------------------|
|  | Consult instructions for use |
|  | Use by Date |