

# Enterovirus/D68 Panel

Enterovirus: RT-qPCR



**Multiplex**

Enterovirus and  
Enterovirus D68

**Optical Channels**

FAM and Cy5

**Fast**

Less than 2 hours  
testing time

The Enterovirus/D68 Panel is an RT-qPCR assay designed for qualitative multiplex in vitro detection of *enterovirus* and *enterovirus D68*.

Positive controls for each target are included in the kit at approximately 10,000 copies/ $\mu$ L.

Optimized Amplification reagents available, see table below for suggested SKU.

For Research Use Only. Not intended for *in vitro* diagnostics. Not approved by the FDA.

**Made in America**

Targets	Optical Channels
Enterovirus	FAM
Enterovirus D68	CY5

**Sample:**

Nasopharyngeal swab or sputum samples.

**Compatibility:**

Bio-Rad CFX96/384, Opus  
Roche LightCycler 480 System  
ABI QuantStudio 3, 5, 6, 7, 12, and 7500 series

## Ordering

Product Name	SKU#	Reactions	Storage	Price
<b>Enterovirus Panel with OneStep Amplification Kit</b>	<b>500020 &amp; 500028</b>	<b>200</b>	<b>-20°C</b>	<b>\$1,150.00</b>
Enterovirus/D68 Panel	500020	200	-20°C	\$200.00
SeqOnce OneStep RT-qPCR Multiplex Amplification Kit 2400	500028	200	-20°C	\$1,100.00